new

ZFT400
Dr. M. Conti

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The root elevator we are introducing is the result of multiyear experience accrued at the maxillofacial surgery unit of the University of Florence during thousands and thousands of advanced surgical extractions.

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
He holds a degree in Medicine and Surgery, perfected in oral surgery, with a specialisation in Dental Surgery. For more than 27 years he has been working at the Maxillofacial Surgery Unit of Florence where he has mainly been involved in outpatient oral surgery. In the field of extractive surgery alone he has carried out more than 5,800 operations. He taught at the specialisation school in oral surgery of the University of Florence, holding the chair of Endodontic Surgery: he has been a lecturer in the annual advanced master in Oral Surgery and a reference tutor for a great number of students. In addition to various scientific papers and participation as a speaker at conventions and congresses, he has recently published an illustrated text on octave surgery for the publishing house, Masterbooks of Florence.

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If we observe the section of the working tip of the majority of root elevators found on the market we can see how they have a "circular arc" shape designed to follow the contours of the root of the tooth, taking it for granted that there is sufficient periodontal space for introducing the elevator.

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
*The use of this shape assumes that there is sufficient periodontal space for being able to insert the elevator. While this is true for simple extractions it is not the case for more complex extractions where, due to the anatomy it is often distorted by the osteotomy and odontotomy. At times the fulcrum searched for is not even in the bone, but instead, inside the tooth itself, as in the case of an odontotomy fissure created intentionally in order to divide the tooth into two or more sections. It is inside these spaces that the active part of the elevator must be engaged and used with rotational movements on its own axis, and only very rarely used as a first-class elevator (see extract of text by Conti Massimo "Chirurgia degli Ottavi" p. 157 – 164 Ed. Masterbooks 2015)**

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The shape of the section of the working tip that we are now proposing is designed to wedge itself inside better and once it is engaged it must be used with small rotational movements along the main axis of the instrument. The flattened back of the active part reduces the height (h) thus shifting the point of the fulcrum (F) to the end and in this way enabling readier and more efficient mobilising compared to the semicircular back that instead has the fulcrum in the middle. Moreover, the moveable fins of the ventral part increase the grip on the tooth.

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The relationship between the length of the handle and the length of the stem makes it possible to avoid areas of shadow created by the operator's hand, also for deep bone cavities induced surgically. These are the considerations that have motivated us to make this new elevator available to dentists in the hope that it will be a helpful aid in their daily practice.