CAD
SYSTEM

LASCOD

Ergominix

A+B

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APPLICATIONS

Control masks for the design of prosthetic structures, creation of provisional devices, reproduction of removable gingival surfaces (XILGUM). Moulding technology using ceramics and composite resins, filling masks for combined prostheses, bite registration keys scannable with CAD systems, models for the repair of prostheses, duplication of plaster models without the use of isolating agents, duplication of whole prostheses. Repositioning and addition of teeth on prostheses and partial dentures, isolation for resin teeth when using the flask method.

CHARACTERISTICS

- Scannable using optical/laser/tactile scanning systems.
- Suitable for moulding methods.
- High details definition (40 micron).
- Final hardness achieved within short time.
- High resistance to compression.
- Compatibility with self- and heat-curing acrylic resins/heat resistance: over 120°C.
- Can be finished using burs.
- High colour contrast.
- High dimensional stability over time.

TECHNICAL DATA	ERGAMIX A+B 70	ERGAMIX A+B 90
Colour	Red	Blue
Base paste/catalyst ratio	1:1	1:1
Mixing time	30"	30"
Working time (23°C/73°F)	2' 45"	2' 45"
Setting time (23°C/73°F)	10'	10'
Final hardness in degrees Shore A	70	90
Elastic recovery	99.5%	99.5%
Deformation on compression	1.8%	0.8%
Size variation (after 24 hours)	-0.02%	-0.02%
Reproduction of details	40 µm	40 µm

The above times can be influenced by:

- variations in the temperature of the material and/or environment (higher temperatures reduce them while lower temperatures prolong them);
- changes in the speed and/or force employed during mixing
- humidity.

INSTRUCTIONS FOR USE

Take equal parts of the two components from the cans using the respective coloured measuring spoons. Comply with the indicated ratio (1:1) to avoid altering the product's final characteristics and the respective working and setting times. Close the cans again as soon as the material has been taken, and do not interchange their respective measuring spoons and caps. Mix vigorously until a uniformly coloured mass without streaks is obtained. Model the mass as needed, adapting it to the relevant surface. Once it has hardened, the product can be removed to proceed with the next phases.

WARNINGS

Wear gloves and protective clothing. **Use only vinyl/nitrile gloves. Avoid direct or indirect contact with latex rubber gloves or with components of such material. In case you have used latex gloves, wash your hands thoroughly. Contact with latex inhibits the hardening of polyvinylsiloxanes.** Store in a cool and dry place. Use within the expiry date and for the applications described.

OTHER AVAILABLE SILICONES

For dental clinics: SILAXIL condensation silicone; KROMOPAN SIL, GHENESYL and OKLUREST addition silicones.

For laboratories: ERGASIL condensation silicone; XILGUM addition silicone.

GUARANTEE

LASCOD SpA guarantees the quality of its products and undertakes to replace them whenever defects in production are found within the expiry dates. LASCOD SpA will not be held responsible and will not recognize any damage that is direct, indirect, accidental or caused by incorrect use of the material. Check the suitability of the product prior to use.

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