CAD  
SYSTEM

**LASCOD**

**Ergominix**

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## APPLICATIONS

Control masks for the design of prosthetic structures, creation of provisional devices, reproduction of removable gingival surfaces (XILGUM). Moulding technology using ceramics and composite resins, filling masks for combined prostheses, bite registration keys scannable with CAD systems, models for the repair of prostheses, duplication of plaster models without the use of isolating agents, duplication of whole prostheses. Repositioning and addition of teeth on prostheses and partial dentures, isolation for resin teeth when using the flask method.

## ADVANTAGES

- Easy and clean mixing thanks to the 1:1 ratio
- Scannable using optical/laser/tactile scanning systems
- Suitable for moulding methods
- High details definition (40 micron)
- Final hardness achieved within short time
- High resistance to compression
- Compatibility with self- and heat-curing acrylic resins/heat resistance: over 120°C
- Can be finished using burs
- High colour contrast
- High dimensional stability over time

## METHOD FOR USE

Take equal parts of the two components from the cans using the respective coloured measuring spoons. Mix vigorously until a uniformly coloured mass without streaks is obtained. Model the mass as needed, adapting it to the relevant surface. Once it has hardened, the product can be removed to proceed with the next phases.

TECHNICAL DATA	ERGAMIX A+B 70	ERGAMIX A+B 90
Colour	Red	Blue
Base paste/catalyst ratio	1:1	1:1
Mixing time	30"	30"
Working time (23°C/73°F)	2' 45"	2' 45"
Setting time (23°C/73°F)	10'	10'
Final hardness in degrees Shore A	70	90
Elastic recovery	99.5%	99.5%
Deformation on compression	1.8%	0.8%
Size variation (after 24 hours)	-0.02%	-0.02%
Reproduction of details	40 µm	40 µm

## WARNINGS/PRECAUTIONS

The times indicated above may be influenced by the following factors:

- Changes to the temperature of the material, the environment and/or the operator's hands during mixing (higher temperatures reduce the times, lower temperatures increase them).
- Changes to the speed and/or force of mixing.
- Humidity.

Comply with the indicated ratio (1:1) to avoid altering the product's final characteristics and the respective working and setting times. LATEX GLOVES INHIBIT THE HARDENING OF ADDITION SILICONES, hence the use of vinyl or nitrile gloves and accurate washing of the hands are recommended. Close the cans again as soon as the material has been taken, and do not interchange their respective measuring spoons and caps. Perfectly clean, degrease and dry the surfaces with which the silicone comes into contact.

## OTHER AVAILABLE SILICONES

**For dental clinics:** SILAXIL condensation silicone; OKLUREST and GHENESYL addition silicones.

**For laboratories:** ERGASIL condensation silicone; XILGUM addition silicone.

## GUARANTEE

LASCOD SpA guarantees the quality of its products and undertakes to replace them if production defects are discovered within the expiry date. LASCOD SpA is not liable and does not admit any claims for any direct, indirect or accidental damages or damages caused by improper use of the material. Check the product's suitability prior to use.

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