
# chromatic alginate

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
Tipo 3: Kromotypo 3; Singletypo 3. Tipo 4: Kromotypo 4; Singletypo 4; Ortotypo.

# chromatic alginate

Measurements taken at 23 °C with deionised water.

## TYPE I: FAST SETTING - TYPE II: NORMAL SETTING

"iQ" does not require special treatment before use. Dust-free material. However, caution is advised when opening the bag and taking the powder. While tap water may be used, the declared times were measured by mixing the powder with deionised water at 23°C.


**Vigorously mix the powder with water.** When using specific mechanical mixers for alginate-based impression materials, follow the manufacturer's instructions.

**Follow the chromatic phase indications which ensure that the product is mixed perfectly and inserted into mouth at the right moment. Violet indicates the mixing time with a spatula, lilac indicates the time available for loading the impression tray, green indicates the time for inserting in the mouth.**

Variations in temperature and/or water hardness and the speed and/or friction during mixing cause small changes in the working time and in the initial setting time, which are automatically indicated by the colour changes. Wash, rinse and shake the impression and pour the model. Keep the container tightly closed in a cool, dry place. Excellent results can be obtained with LASCOD plasters:

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