
Kromoglass

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Water-mixable glass ionomer cement for luting.

LASCOD Spa - Via L. Longo, 18 - 50019 Sesto F.no (Firenze), Italy
 tel. (+39) 055 421 57 68 (r.a.) - lascod.italy@lascod.it
 www.lascod.com

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KROMOGLASS 3

Water mixable glass ionomer cement for luting. Compliant with ISO 9917-1 specification. For professional use.

USES: permanent cementing of bridges, crowns, inlays, onlays and orthodontic bands. Cementation of ceramic restorations (silicate ceramics, zirconium oxide, aluminum oxide).

PREPARATION OF THE AREA: Kromoglass 3 is an extra fine glass ionomer cement with extraordinary adhesion qualities to enamel and dentine. For best results, the enamel and dentine must be cleaned and not excessively dried to prevent postoperative sensitivity. Thoroughly dry the surfaces of the restoration.

PREPARATION AND USE OF THE CEMENT: introduce water, preferably distilled, into the dropping bottle by aspiration. Shake the powder before use. For each level measure of powder use one drop of water. After taking the powder close the container tightly with its own cap. Remove any remaining product from the measuring spoon after use. Hold the dropper in a vertical position for dosing the water. Pour the powder and the water onto a glass slide and, using an abrasion resistant spatula made of plastic or metal, mix the powder with water starting off with small quantities and then gradually increase the amounts. Carry out the preparation at a temperature of 23 +/-1°C. With different temperatures use a glass slide with a temperature included between those mentioned above. Higher temperatures will shorten working time, while lower temperatures will lengthen it. The right consistency is obtained when the prepared mixture can flow but still sticks to the spatula. Fit the restoration into place taking care to exert the right amount of pressure and not to use too much cement. During the application and the setting process make sure that the working area is kept dry. Avoid using dentine insulating materials in order to maintain the fixing capacity of the cement intact. Clean the instruments used in water before the adhesive cement sets.

Mixing ratio	1.0 g of water to 3.3 - 3.8 g of powder
Mixing time*	30" - 40"
Working time, including mixing time*	2'30" - 3'30"
Net setting time**	3'30" - 5'

* Measurements taken at 23 +/-1°C. ** Measurements taken at 37 +/-1°C.

WARNING: For use, wear protective gloves, eye protection and protective clothing. Use within the expiry date and for the applications described. Do not use in case of hypersensitivity to: fluorosilicate glass, polyacrylic acid, tartaric acid. As generally recognized for glass ionomer cements, the product is slightly cytotoxic. Store in a cool and dry place. The product fears humidity. Do not use in case of aggregation of the powder.

OTHER GLASS IONOMER CEMENTS: Kromoglass 2 (filling cement).

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