Kromoglass 2

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Water-mixable glass ionomer cement for fillings. Radiopaque.

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KROMOGLASS 2
 Water mixable glass ionomer cement for fillings. Radiopaque.
 Compliant with ISO 9917-1 specification. For professional use. (VITA shade scale - A3)

USES: cavity restoration in deciduous teeth especially of class 1. Repair of cavities with carious lesions, repair of erosions of the enamel and the roots in the neck area of the tooth. Permanent Class 3 fillings, inclusive Class 5 fillings and small Class 1 fillings. Fillings of enlarged fissures and supporting fillings for crowns, bridges and base fillers.

PREPARATION OF THE CAVITY: the creation of constraints is not normally necessary due to the excellent adhesion of this cement to both the dentine and the enamel, with which it creates chemical bonds. Undercuts can nevertheless give the filling additional hold. The thickness of the filling should not be less than 1 mm. Where there are carious lesions, the cavity must be cleaned with 3% hydrogen peroxide, thoroughly rinsed with water and carefully dried. If the thickness of the residual dentine exceeds 1.5 mm additional protection with calcium hydroxide will not be necessary. In order not to impair the adhesive ability of the cement, avoid the use of an insulating material when not strictly necessary.

PREPARATION AND USE OF THE CEMENT: introduce water, preferably distilled, into the dropping bottle by aspiration. Shake the powder before use. For each level measure of powder use one drop of water. After taking the powder close the container tightly with its own cap. Remove any remaining product from the measuring spoon after use. Hold the dropper in a vertical position for dosing the water. Pour the powder and the water onto a glass slide and, using an abrasion resistant spatula made of plastic or metal, mix the powder with water starting off with small quantities and then gradually increase the amounts. Carry out the preparation at a temperature of 23+/-1°C. With different temperatures use a glass slide with a temperature included between those mentioned above. Higher temperatures will shorten working time, while lower temperatures will lengthen it. The right consistency is obtained when the prepared mixture can flow but still sticks to the spatula. Avoid bevelling the edges. During the application and the setting process make sure that the working area is kept dry. The use of rubber dams is advisable. After the application of Kromoglass 2 immediately cover the filling with a mould, as it is sensitive to humidity at this stage. After removing the mould, the filling must be covered with a suitable varnish. 5-6 minutes after the completion of setting, excess cement can be removed and the levels checked. The best results of the final cleaning (flexible discs or fine grain diamond) are obtained after 24 hours. Clean the instruments used in water before the adhesive cement sets.

Mixing ratio	1.0 g of water to 4.8 - 5.2 g of powder
Mixing time*	30" - 40"
Working time, including mixing time*	2'30" - 3'30"
Net setting time**	3' - 4'30"

* Measurements taken at 23 +/-1°C. - ** Measurements taken at 37 +/-1°C.

WARNING: The product can cause eye injury. For use, wear protective gloves, eye protection and protective clothing. Use within the expiry date and for the applications described. Do not use in case of hypersensitivity to: fluorosilicate glass, polyacrylic acid, tartaric acid. As generally recognized for glass ionomer cements, the product is slightly cytotoxic and, if used on low dentin thicknesses without adequate protection, it can have effects on the pulp (irritation and /or necrosis). Store in a cool and dry place. The product fears humidity. Do not use in case of agglomeration of the powder.

OTHER GLASS IONOMER CEMENTS: Kromoglass 3 (luting cement).

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