CAD
SYSTEM

LASCOD

Xilgum

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ADVANTAGES

- Easy and safe handling
- Maximum fluidity
- Rapid setting time
- High dimensional stability over time
- Suitable for direct and indirect methods
- No retraction / deformation
- Scannable using optical/laser/tactile scanning systems
- Easy to finish with cutters or burs

USE - INDIRECT METHOD - 1. Before cutting the model, prepare a mask with a thickness suitable for covering the involved part. - 2. Remove the gingival area to be reproduced from the plaster model (a thickness of at least 2 mm is recommended). - 3. Use a burr to cut at least two holes (inlet/outlet) into the silicone mask. - 4. Accurately clean the plaster model. Internally isolate the silicone mask via the application of SEP FLUID. Reposition the mask on the model. - 5. Slowly inject XILGUM starting on the vestibular side, applying a constant pressure on the dispenser's lever, until the material comes out of the hole on the opposite side of the mask. To prevent the formation of bubbles, do not stop the extrusion of XILGUM until the mask has been completely filled. Wait for the material to harden completely. - 6. Remove the silicone mask and the reproduced gingiva. During this operation we recommend to pull out or cut the injection channels to prevent tearing the area reproduced with XILGUM. - 7. Proceed with finishing. To eliminate excess material use cutters or tungsten burs.

USE - DIRECT METHOD - 1. Clean and dry the impression. - 2. Uniformly apply the isolating liquid SEP FLUID onto the entire involved surface and wait for it to dry. - 3. Directly extrude the material onto the relevant area by applying a constant pressure on the dispenser's lever. It is recommended to keep the tip immersed in the material to prevent the inclusion of air bubbles. The material's fluidity ensures a homogenous distribution even in areas that are difficult to reach. - 4. Wait for the material to harden completely. - 5. Create the model in your normal manner using your usual materials. It is not necessary to isolate the exposed XILGUM surfaces. - 6. Once the material used has set, separate the artefact from the impression by pulling out the gingiva reproduced with XILGUM. - 7. Finish using cutters or burs.

TECHNICAL DATA	XILGUM	XILGUM SOFT
Base paste/catalyst ratio	1:1	1:1
Working time (23°C/73°F)	2'	1' 45"
Setting time (23°C/73°F)	10'	6'
Reproduction of details	20 µm	20 µm
Shore A hardness	72	40
Size variation (after 24 hours)	-0,02 %	-0,02 %
Elastic recovery	99,9 %	99,9 %

WARNINGS/PRECAUTIONS - The times indicated above may be influenced by the following factors:

- Changes to the temperature of the material and of the environment (higher temperatures reduce the times, lower temperatures increase them).
- Humidity.

Perfectly clean, degrease and dry the surfaces with which XILGUM comes into contact. For correct storage of the remaining product in the cartridge leave the mixing tip installed until the next use. Before installing a new mixing tip, make sure that the silicone properly flows from the nozzles.

XILGUM SEP FLUID - SEPARATOR LIQUID

XILGUM SEP FLUID is a ready-to-use separator agent that permits isolation between the surfaces of materials of the same nature to prevent adhesion. It is used to create gingival reproductions on models in association with XILGUM and more generally in all cases where adhesion between compatible materials must be prevented. XILGUM SEP FLUID may be applied to addition silicones (ERGAMIX 70-90 Shore A) and condensation silicones (ERGASIL). Only with polyether-based materials is it necessary to wait at least 5 hours after the impression has set.

INSTRUCTIONS FOR USE - The product may be applied directly onto the dry impression (direct method) or onto silicone masks (indirect method). Make sure that the surfaces to be isolated are perfectly clean and dry. Uniformly spread XILGUM SEP FLUID with the brush attached to the cap and wait for it to dry completely. Drying may be accelerated by a gentle air stream.

WARNINGS - Close the bottle well after use.

OTHER AVAILABLE SILICONES - For dental clinics: SILAXIL condensation silicone; KROMOPAN SIL, GHENESYL and OKLUREST addition silicones. For laboratories: ERGASIL condensation silicone; ERGAMIX addition silicone.

GUARANTEE - LASCOD SpA guarantees the quality of its products and undertakes to replace them if production defects are discovered within the expiry date. LASCOD SpA is not liable and does not admit any claims for any direct, indirect or accidental damages or damages caused by improper use of the material. Check the product's suitability prior to use.

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