

ITALIAN
QUALITY
ALGINATE

Value for Money



Now with the
**NEW
MANGO
FLAVOUR!**



LASCOD

iQ will suit your expectation for Value at the right Price!



- Versatile and easy to use with the chromatic indicator (**chrom and tropical**).
- You can choose! Fast or Normal Setting.
- Thixotropic
- High definition details (**50µm**).
- Extremely elastic (**recovery from deformation 95.5%**).
- Remember! Minimum **60" in mouth!**
- Easy to read final colors.
- Compatible with disinfectants liquids.
- Long expiry date (**5 years**)
- No lead or heavy metals content.
- Dust free



Type 1	IQN302	IQM302	IQO302
Type 2	IQU302	IQT302	IQV302
Working time* - Type 1		1' 50"	
Working time* - Type 2		2' 50"	
Minimum time in mouth (37°)		60"	
Accuracy (µm)		50	
Elastic recovery (%)		95.5	
Flavour	mint	mint	mango
Colour	green		

* Measurement at 23°C - Mixing time included.



iQ box
contains
24x450g. bags
(10,8 kg)

The photographic images here reproduced are purely indicative and are not necessarily identical to the actual products.

LASCOD

For further information activate the QR Code reader on your smartphone.

